## DOCUMENTO

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1. Mann J, Truswell AS. Effect of controlled breakfast on serum cholesterol and triglycerides. *Am J Clin Nutr* 1971; 24: 1300.
2. Peric-Golia L. Postprandial serum-cholesterol. *Lancet* 1972; ii: 876.
3. Wood PD, Haskell WL. The effect of exercise on plasma high-density lipoproteins. *Lipids* 1979; 14: 417-427.
4. Fuller JH, Grainger SL, Jarret RJ, Keen H. Possible seasonal variation of plasma lipids in a healthy population. *Clin Chim Acta* 1974; 52: 305-310.
5. Warnick GR, Albers JJ. Physiological and analytical variation in cholesterol and triglycerides. *Lipids* 1976; 11: 203-208.
6. Van Gent CM, Van Der Voot H, Hessel LW. High-density lipoprotein cholesterol, monthly variation and association with cardiovascular risk factors in 1000 forty year-old Dutch citizens. *Clin Chim Acta* 1978; 88: 152-176.
7. Glomset JA, Norum KR. The metabolic role of lecithin: cholesterol acyl-transferase: Perspectives from pathology. *Adv Lipid Res* 1973; 11: 1-65.
8. Zilversmit DB, Hughes LB, Balmer J. Stimulation of cholesterol ester exchange by lipoprotein-free rabbit plasma. *Biochim Biophys Acta* 1975; 409: 393-398.
9. Chajek T, Fielding CJ. Isolation and characterization of human serum cholesteryl ester transfer protein. *Proc Natl Acad Sci USA* 1978; 75: 3445-3449.
10. Bachorik PS, Walker R, Brownell KD, Stunkard AJ, Kwiterovich PO. Determination of HDL-cholesterol levels in stored human plasma samples. *J Lipid Res* 1980; 21: 608-616.
11. Wood PD, Bachorik PS, Albers JJ, Stewart CC, Winn C, Lippel K. An investigation of the effects of sample aging on total cholesterol values determined by the automated ferric chloride-sulfuric acid and Lieberman-Burchard procedures. *Clin Chem* 1980; 26: 592-597.
12. Lewis B. Lipids, en «Chemical diagnosis of disease» Brown SS, Mitchell FL, Young DS (editores) 1979; Elsevier/North-Holland pp 287.
13. Tan MH, Wilmhurst EG, Gleason RE, Soeldner JS. Effect of posture on serum lipids. *N Eng J Med* 1973; 289: 416-419.
14. Statland BE, Bokelund H, Winkel P. Factors contributing to intra-individual of serum constituents: 4. Effects of posture and tourniquet application on variation of serum constituents in healthy subjects. 1974; 20: 1513-1519.
15. Koerselman HB, Lewis B, Pilkington TRE. The effect of venous occlusion on the level of serum cholesterol. *J Atheroscler Res* 1961; 1: 85-88.
16. Page JH, Moinuddin M. The effect of venous occlusion on serum cholesterol and total protein concentration- a warning. *Circulation* 1962; 25: 651-652.
17. Nishida J, Kummerow FA. Interaction of serum lipoproteins with hydroperoxide of methyl linoleate. *J Lipid Res* 1960; 1: 450-458.
18. Nichols A, Rehnberg CS, Lindgren FT. Gas chromatographic analysis of fatty acids from dialyzed lipoproteins. *J Lipid Res* 1961; 2: 203-207.
19. Hatch FT, Lees RS. Practical methods for plasma lipoprotein analysis. *Adv Lipid Res* 1968; 6: 1-68.
20. Schuh J, Faircloth GF, Haschemeyer RH. Oxygen mediated heterogeneity apo-low density lipoprotein. *Proc Natl Acad Sci Usa* 1978; 75: 3173-3177.
21. Schwertner HA, Friedman HS. Changes in lipid values and lipoprotein patterns of serum samples contaminated with bacteria. *Am J Clin Pathol* 1973; 59: 829-835.
22. Knoetgen V, Petek W, Kostner G. Alteration of serum lipid values and lipoproteins by bacteria as a possible cause of artifacts in the screening of dyslipoproteinemias. *Klin Wochenschr* 1976; 54: 643-645.
23. Krishnaiah KV, Wiegandt H. Demonstration of a protease like activity in human serum low density lipoprotein. *FEBS Lett* 1974; 40: 265-268.
24. Chapman MJ, Kane JP. Stability of the apoprotein of human serum low density lipoprotein: Absence of endogenous endopeptidase activity. *Biochem Biophys Res Commun* 1975; 66: 1030-1036.
25. Bachorik PS, Albers JJ, Ellefson RD, Kane JP. Collection of blood samples for lipoprotein analysis. *Clin Chem* 1982; 28: 1375-1378.
26. Del Gatto L, Nichols AV, Lindgren FT. Effects of freeze-thawing and storage on ultracentrifugal properties of human serum lipoproteins. *Proc Soc Exp Biol Med* 1959; 101: 59-61.
27. Mills GL, Wilkinson PA. Some effects of storage on plasma lipoproteins. *Clin Chim Acta* 1962; 7: 685-693.
28. Warnick GR, Albers JJ. High density lipoprotein cholesterol (HDL CH) quantitation: Effect of plasma storage on heparin-Mn<sup>2+</sup> supernatant. *Clin Chem* 1979; 25: 1908.
29. Warnick GR, Mayfield C, Albers JJ. Evaluation of quality control materials for high density-lipoprotein cholesterol quantitation. *Clin Chem* 1981; 27: 116-123.